SIMPLE Ag SARS-CoV-2 / Flu A+B

USO EXCLUSIVO PROFESIONAL

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SIMPLE Ag SARS-CoV-2 / Flu A+B

FOR PROFESSIONAL USE ONLY

PLEASE READ THESE INSTRUCTIONS CAREFULLY BEFORE USING THE TEST.

INTRODUCTION

• Intended use

A single-step qualitative immunochromatographic test for simultaneously detecting SARS-CoV-2, Influenza A, and Influenza B virus antigens in nasopharyngeal and oropharyngeal swab type human respiratory samples using separate strips.

• General information

The SARS-CoV-2 virus is a new coronavirus that causes an infection known as COVID-19. Transmission is mainly via aerosols and drops emitted when an infected person speaks, breathes, coughs, or sneezes. This is a very contagious virus, so rapidly diagnosing the disease is particularly important with a view to isolating the patient and preventing the virus from propagating.

Influenza A and B virus infection is also mainly a respiratory condition, coursing with similar symptoms to COVID-19, thereby supporting their separate simultaneous detection using this test.

The Ag SARS-CoV-2 / Flu A+B test can detect both the SARS-CoV-2 virus nucleoprotein, and the Influenza A and B virus nucleoproteins in the same sample.

This test is exclusively for *in vitro* diagnostic use.

The results of this test are intended to support data derived from a clinical evaluation of the patient; the test does not constitute a definitive diagnosis in itself. A negative result does not rule out an infection by SARS-CoV-2 or Influenza. In the event of obtaining a negative result for a patient presenting the typical symptomatology from these types of infections,

confirming the result using RT-PCR is recommended. Any decisions (treatment, isolation, etc.) should not be based solely on the test. This test is for professional use. This is not a self-diagnostic test.

Intended test population

The test is intended for use on the general population, as any individual is susceptible to infection by these two viruses.

Disease or infection incidence among the intended test population

SARS-CoV-2: The disease incidence among the worldwide population is constantly changing as it depends on how the pandemic progresses and the region being analysed at any given time. Two years from the onset of the pandemic, there have been more than 500 million cases around the world.¹

Influenza: worldwide, from the start of November 2019 to the end of December 2020, 614,907 respiratory samples collected during hospital and outpatient surveillance were positive for Influenza (19% of all analysed samples). Of these positive samples, 63% and 37% were typified as influenza A and influenza B respectively.² The WHO estimates that annual flu epidemics cause approximately 3 to 5 million cases of serious disease, and 300,000 to 500,000 deaths.

Characteristics of the virus and its infection

SARS-CoV-2: coronaviruses are single-stranded RNA viruses that belong to the Coronaviridae family. They have an external envelope that looks like a crown (corona), hence the name. The four most important structural proteins are the spike protein (Spike - S-), nucleocapsid protein (N), membrane protein (M), and envelope protein (E). Of these, the S protein is responsible for interacting with human cell receptors, making it an important target for neutralising antibodies and vaccine development.³

Coronaviruses are classified into three types: α , β (viruses only pathogenic to mammals) and γ (mainly infecting birds). Seven types of human coronaviruses have been identified to date: HCoV-229E, HCoV-NL63, HCoV-OC43, HCoV-HKU1, SARS-CoV, MERS-CoV, and the new coronavirus first identified in Wuhan in 2019, SARS-CoV-2, which causes the COVID-19 disease. The new coronavirus is associated with cases of pneumonia with common symptoms including dry cough, fever, fatigue, or loss of smell. It can easily progress to severe pneumonia, acute respiratory distress syndrome, multiple organ failure, and even death.^{4,5} It has a mortality rate of 1-2%, depending on the geographical location.

Influenza: the Influenza virus is an RNA virus containing the hemagglutinin (H) and neuraminidase (N) proteins in its envelope. The H protein is responsible for interaction with host cells, and therefore, the virus' entry, while the N protein participates in both virus replication and its exit from the cell. The H and N viral proteins show most variation in terms of antigens, while the virus'

nucleoprotein, detected by the antibodies used for both Influenza A and Influenza B, is much more conserved. Influenza A viruses are divided into subtypes based on the H and N proteins on the virus' surface. The influenza A virus subtypes currently circulating more commonly among people are: A (H1N1) and A (H3N2). Influenza B viruses, on the other hand, are classified into two lineages: B/Yamagata and B/Victoria.

Symptoms associated with flu virus infection vary from a mild respiratory disease - limited to the upper respiratory tract and characterised by fever, throat pain, nasal secretion, cough, headache, muscle pain and fatigue - to severe pneumonia, which in some cases can prove fatal, whether due to the influenza virus itself or a secondary bacterial infection in the lower respiratory tract. In some cases, infection by the Influenza virus can also cause a wide range of non-respiratory complications that affect the heart, central nervous system, and other organ systems.⁶

THE BASIS OR BASIC PRINCIPLES OF THE TEST

The Ag SARS-CoV-2/Flu A+B test is an immunochromatographic test with two strips inside a double housing.

1. **The Ag SARS-CoV-2** strip uses a combination of:

a. Purple latex particles conjugated to an antigen recognised by a specific antibody for said antigen bound to the membrane, comprising the so-called test control line (position C on the housing).

b. Red latex particles conjugated to a specific antibody that targets the SARS-CoV-2 virus nucleoprotein acting in conjunction with another specific antibody for the protein located on the test's membrane, below the control line (position T on the housing).

2. **The Influenza strip** uses a combination of:

a. Purple latex particles conjugated to an antigen recognised by a specific antibody for said antigen bound to the membrane, comprising the so-called test control line (position C on the housing).

b. Blue latex particles conjugated to a specific antibody that targets the Influenza type B virus nucleoprotein acting in conjunction with another specific antibody for said protein located on the membrane, below the control line (position T2 on the housing).

c. Red latex particles conjugated to a specific antibody that targets the Influenza type A virus nucleoprotein acting in conjunction with another specific antibody for said protein located on the membrane, below the line that detects Influenza type B (position T1 on the housing).

This test first involves treating the sample with the sample dilution buffer (included in the kit) to extract the SARS-CoV-2 and Influenza viruses. Once extracted, simply add a specific volume of the extract to the reagent strips and wait for 15 minutes.

The coloured particles migrate when the extracted sample flows through the strips. In the event of a

positive sample, the antibodies present in the particles will react with the virus antigens, forming a particle-antigen complex, which will also be retained on the membrane by other specific antibodies for that antigen. Different coloured lines will appear depending on the virus content in the sample. These lines are used to interpret the result after 15 minutes of incubation at room temperature (see Fig. 1).

These instructions for use apply to the trade references for the product: 9.242.XXX.YY.ZZZ.

MATERIALS INCLUDED IN THE KIT

- 20 reaction devices.
 - 20 sterile swabs
- Goodwood ref. GW-1237NP  0197
- 1 bag of 20 individual vials containing sample dilution buffer.
 - Positive controls and the corresponding instructions for use in the relevant trade reference.
 - Data matrix in the relevant trade reference.
 - Instructions for use.

The product's safety sheet contains information about the composition of the reagents. An electronic copy of the sheet can be requested via email to msds@operon.es

MATERIALS NOT INCLUDED IN THE KIT

- Stopwatch

Controls are available for the user to validate the results, and can also be purchased as a separate trade reference.

EQUIPMENT FOR INTERPRETING THE RESULTS

The test can be interpreted visually or using OPERON's immunochromatographic test reader.

For this diagnosis, there is evidence of a 99% agreement between visual interpretation and the OPERON IC strip reader.

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PRECAUTIONS

Check all the applicable precautions from this example list:

1. Patient samples must be handled with care as they may contain infectious agents. Disposable gloves and other protective measures that may eventually be required must be used throughout the process.
2. The sample dilution buffer contains sodium azide as an anti-microbial agent. Avoid direct contact with skin and mucous membranes. Dispose of appropriately. Do not use the buffer if there are any signs of contamination or precipitation.

3. Do not store or prepare food, eat, drink, or smoke in the area where the reagents and samples are handled.

4. On completing the task, clean all work surfaces with soap and water, and disinfect them with a suitable solution. Lastly, dispose of the protective gloves correctly and wash your hands with soap and water, rubbing them well.

5. Do not exchange components from kits with different batch numbers.

6. If the test is stored refrigerated, allow all the kit's components and the samples to reach room temperature, as cold reagents and/or samples can reduce the test's functionality. Allowing 20 to 30 minutes to reach room temperature is recommended.

7. Only use the reagents *in vitro*.

8. Do not use any of the kit's components after the use-by-dates.

9. If the outer box is damaged, the product can still be used providing none of the components have been damaged.

10. Throw the test away if the primary packaging is either broken or impaired.

11. Remove the strip or device from the aluminium bag when the test is going to be performed to avoid unnecessary excessive exposure of the product to any potentially damaging environmental factors.

12. Adding the correct extracted sample volume to both sample addition zones indicated on the housing is very important. If lower than indicated, the chromatography may not process correctly because insufficient sample reaches the reaction zone. If higher than indicated, the chromatography may not process correctly due to the sample volume surpassing the strip's absorption capacity.

13. Please dispose of the used product in accordance with current legislation.

14. Do not use the test if there are any coloured lines in the results zone prior to use.

15. Mix all respiratory samples thoroughly before performing the test in order to ensure a representative sample is used.

16. Avoid sampling from any viscous zones remaining in the sample after mixing. If the pipette blocks, return the sample to the container and sample it again.

17. Do not use sputum or saliva samples as invalid results may be obtained.

19. Do not discard the kit's outer box until all the content has been used. The box contains essential information about the product's CE marking and the batch number.

20. Any serious incident involving the product must be reported to the manufacturer and the relevant authority in the member state where the user and/or patient are located.

STORAGE

The Simple Ag SARS-CoV-2/FluA+B product can be

stored at any temperature from +2 to +30 °C. The use-by date of each component is printed on the aluminium wrappers. The reagent strips must be used once the protective packaging has been opened after reaching room temperature (if stored refrigerated). The product maintains the claimed stability once open.

SAMPLES

- The test is designed to analyse nasopharyngeal and oropharyngeal swab type human respiratory samples.
- Samples suspended in the product's dilution buffer have to be analysed immediately after sampling, storing them for later use is not recommended.
- Processing the samples as soon as possible after collection is recommended. For samples collected in a different viral transport medium to the product's sample dilution buffer, please see the specifications for each medium when it comes to the stability of suspended samples and storage conditions.
- Limit sample freeze/thaw cycles, as they could negatively affect the results.
- Do not centrifuge samples prior to using them with the Simple Ag SARS-CoV-2/Flu A+B test, as eliminating cell material could adversely affect the test's sensitivity.
- Several swab types for sampling have been evaluated. If the sample is nasopharyngeal, using the Goodwood brand nasopharyngeal swab included in the kit is recommended, or the Copan Flock Technologies FLOQSwabs (ref. 503CS01) swab with similar characteristics. For oropharyngeal samples, using Copan Flock Technologies FLOQSwabs oropharyngeal swabs (ref. 502CS01) is recommended.
- The test is compatible with at least the Vircell and CPM Scientifica brand transport mediums. If using different transport media to the product's dilution buffer, the test's characteristics may vary depending on the type and volume used.

PREPARING SAMPLES AND THE TEST PROCEDURE

General Note: use disposable gloves throughout the test as infectious samples are being handled. On completion, follow the hygiene procedures in point 4 of the "Precautions" section.

THE TYPE OF SAMPLE TO USE DETERMINES THE PROCEDURE TO FOLLOW

- **Procedure for using the test with oropharyngeal/nasopharyngeal swab samples in the test dilution buffer.**
 1. For nasopharyngeal samples, use the swab included in the kit to collect the sample. For oropharyngeal samples, using the swabs indicated in the "SAMPLES" section is recommended, or others with similar characteristics should they not be available.
 2. Introduce the swab soaked in sample into the vial containing the sample dilution buffer supplied with the kit.

3. Vigorously twist the swab in the buffer for 30 seconds to extract the sample correctly.
4. Remove the swab, pressing it against the tube walls to ensure as much sample as possible is released.
5. Dispose of the swab.
6. Put the dropper cap into the tube and wait for 2 minutes.
7. Remove the device from its packaging and place it on a flat surface.
8. Gently shake the vial-dropper containing the extracted sample without tipping it over.
9. Slowly add 5 drops to the sample addition zone on each strip (rectangular windows marked with an arrow), allowing each drop of sample to absorb first.
10. Wait **15 minutes**, then read and interpret the result.

- **Procedure for using the test with oropharyngeal/nasopharyngeal swab samples in viral transport media.**

1. Place 160 μL (or 7 drops) of the sample dilution buffer into a plastic microtube (not included with the kit).
2. Add 160 μL of sample using a micropipette (not included with the kit) and mix by pipetting, ensuring not to create foam (avoid using a vortex mixer).
3. Remove the reaction device from its packaging and place it on a flat surface.
4. Add 140 μL of the sample+dilution buffer mix with a micropipette (not included with the kit) to the strip's sample addition zone (a rectangular window marked with an arrow). Add slowly, drop by drop, to ensure the sample absorbs correctly.
5. Wait **15 minutes**, then read and interpret the result.

READING THE RESULTS


The strips shown in Fig. 1 are an example of the different results that can be obtained with the Ag SARS-CoV-2/Flu A+B test.

Different coloured lines can appear on each reagent strip, bordered by black horizontal lines printed on both sides of the housing.

- **The Ag SARS-CoV-2 strip** (located on the left of the housing).
 - a. Purple line: this is the control line, which should always appear, as it indicates that the chromatography worked correctly.
 - b. Red line: this indicates that the sample contains the SARS-CoV-2 virus.
- **The Influenza strip** (located on the right of the housing).
 - a. Purple line: this is the control line, which should always appear, as it indicates that the chromatography worked correctly.
 - b. Blue line: this indicates that the sample contains influenza type B virus.
 - c. Red line: this indicates that the sample contains influenza type A virus.

The purple control line must always appear.

Figure 1 shows examples of possible results. The positive results shown represent individual viral infections. In the case of coinfections, different coloured lines will be shown, red-blue or red-red, on one or both strips).

Strip 1: **NEGATIVE** result

Only a **PURPLE** horizontal line appears on both strips, aligned with the letter "C" on the housing. This is the control line, which must always appear, as it indicates that the chromatography worked correctly.

Strips 2-4: **POSITIVE** results

Strip 2: **SARS-CoV-2 POSITIVE** - Influenza A and B Negative

On the Ag SARS-CoV-2 strip (located on the left of the housing), a **PURPLE** (control) line appears aligned with the letter "C", along with a **RED** line aligned with the letter "T".

On the influenza strip, only the **PURPLE** control line appears, aligned with the letter "C".

Strip 3: **SARS-CoV-2 Negative** - **Influenza A POSITIVE** - Influenza B negative

On the Influenza strip (located on the right of the housing), a **PURPLE** (control) line appears aligned with the letter "C", along with a **RED** line aligned with the symbol "T1".

On the Ag SARS-CoV-2 strip (located on the left of the housing), only a **PURPLE** (control) line appears aligned with the letter "C".

6. Strip 4: **SARS-CoV-2 negative** - **Influenza A negative** - **Influenza B POSITIVE**

On the Influenza strip (located on the right of the housing), a **PURPLE** (control) line appears aligned with the letter "C", along with a **BLUE** line aligned with the symbol "T2".

On the Ag SARS-CoV-2 strip (located on the left of the

housing), only a **PURPLE** (control) line appears aligned with the letter "C".

Strip 5: **INVALID** result

The control lines do not appear or the purple colour of the control lines is clearly altered. This indicates that the test did not work correctly. In the event of an invalid result, repeating the test using a new reaction device is recommended, strictly adhering to the instructions for use described in this manual.

Any line that appears after 15 minutes of reaction time will have no diagnostic value.

NOTE: a doctor is responsible for establishing the final and definitive diagnosis. This test only detects the SARS-CoV-2 and influenza respiratory viruses in a sample. It does not constitute an argument to state that a person is infected with those viruses.

The test can be interpreted visually or by using OPERON's immunochromatographic test reader. The quick guide for interpreting the test using that reader is available from the following website address: <https://operon.es/es/operon-lateral-flow-reader/>
Important: check that no dust particles, fibres, etc. are present in the results window prior to introducing the test into the reader, as they can interfere with the reading.

QUALITY CONTROL

The test is invalid should no (coloured) lines appear, either because it was performed incorrectly or because the reagents have deteriorated. In this case, repeat the analysis, strictly following the operating protocol detailed in these instructions for use.

WARNING: Including our controls with a known result is recommended to ensure that correct data is obtained.

LIMITATIONS OF THE PROCEDURE

1. The Ag SARS-CoV-2/Flu A+B test can (separately) identify SARS-CoV-2 and Influenza A and B viruses, detecting their presence in human nasopharyngeal/oropharyngeal swab type respiratory samples, as long as the viral load is equal to or above the assay's limit of detection.

2. This test is qualitative, not quantitative, although the intensity of the positive lines do relate to the amount of detectable analyte in the sample.

3. A large amount of respiratory samples were evaluated to ensure the test functions correctly. The results correlated well against other reference techniques (such as RT-PCR Real Time). However, this study does not rule out any potential interference to the test's functionality when analysing other respiratory samples.

4. Monitoring the reaction time is important. If the reaction time is shorter than indicated, samples containing a quantity of analyte above the limit of sensitivity will be easy to see, but those on the limit will not appear. If the reaction time is longer than indicated, the test's sensitivity could be affected and may give rise to incorrect results.

5. The test results must be interpreted alongside the available data from epidemiological studies, the patient's clinical evaluation, and other diagnostic procedures.

6. Adding the correct amount of sample to the reagent strips is very important (5 drops using the plastic pipette supplied with the kit), as a smaller amount may cause false negative results, whereas a larger amount could prevent the chromatogram from developing correctly.

7. Unsuitable collection, transportation, and/or handling of the sample could impair the results.

8. In addition to a very intense blue signal on the test zone corresponding to Flu B, if samples have a very high load of Influenza B, an undefined, fainter blue colour may appear on the test zone for Flu A. This does not mean an Influenza A virus infection is present.

9. The antibodies present in the test may not detect all antigenic variants of the SARS-COV-2 and influenza viruses, or any new strains that may appear over time. The Analytical Sensitivity section lists all the strains analysed with the Ag SARS-COV-2/Flu A+B test.

ANALYTICAL SENSITIVITY

▪ Ag SARS-COV-2 strip

The analytical sensitivity with recombinant N protein is 1.56 ng/mL and the limit of detection with the SARS-CoV-2 (USA-WA1/2020, ATCC ref. NR-52287) virus is:

- 3.1 x 10² TCID₅₀/mL in a nasopharyngeal sample.

- 1.1 x 10² TCID₅₀/mL in an oropharyngeal sample.

With an analytical sensitivity similar to when detecting the originally isolated SARS-CoV-2 virus, the test detects the following variants:

– the Gama P.1 variant (Brazil), Sino Biological, Cat. n. 40588-V07E11, N protein expressing the P080R mutation.

– the Delta B.1.617.2 variant (India), Sino Biological, Cat. n. 40588-V07E2, N protein expressing the D63G, R203M and D377Y mutations.

– the Beta B.1.351 variant (South Africa), Sino Biological, Cat. n. 40588-V07E9, N protein expressing the T205I mutation.

– the Omicron B.1.1.529 variant (South Africa), Cat. n. 40588-V07E34, N protein expressing P13L, R203K, and G204R mutations, and ERS31-33 deletion.

– the Beta B.1.351 variant, SARS-CoV-2 virus USA/MD-HP01542/2021 irradiation inactivated, BEI Resources, Cat. n. NR-55351.

– the Alpha variant, B.1.1.7, SARS-CoV-2 virus hCoV-19/USA/CA_CDC_5574/2020 heat inactivated, BEI

Resources, Cat. n. NR-55245.

▪ INFLUENZA strip

Type A:

- Analytical sensitivity with Influenza type A (H1N1), strain A/New Caledonia/20/99 IVR116 chemically inactivated: 8 ng/mL.

- Analytical sensitivity with Influenza type A (H1N1), strain A/New Caledonia/20/99 heat inactivated: 19.5 TCID₅₀/mL.

- Analytical sensitivity with Influenza type A (H1N1), strain A/Wisconsin/67/2005 (viable virus): 3.6 x 10⁵ CEID₅₀/mL.

Type B:

- Analytical sensitivity with Influenza type B, lineage Yamagata, strain B/Florida/07/04 chemically inactivated: 16 ng/mL.

- Analytical sensitivity with Influenza type B, lineage Yamagata, strain B/Florida/07/04 heat inactivated: 2.5 x 10³ TCID₅₀/mL.

- Analytical sensitivity with Influenza type B, lineage Yamagata, strain B/Florida/07/04 (viable virus): 3.6 x 10⁵ CEID₅₀/mL.

The Influenza strip analytical sensitivity was analysed with other variants of the virus (A and B):

Strain	Subtype/lineage	Analytical sensitivity
A/New Caledonia/20/99	H1N1	8 ng/mL
A/Solomon Island/03/06	H1N1	31 ng/mL
A/Beijing/262/95	H1N1	31 ng/mL
A/Taiwan/1/86	H1N1	16 ng/mL
A/Brisbane/10/07	H3N2	8 ng/mL
A/Kiev/301/94	H3N2	16 ng/mL
A/Panama/2007/99	H3N2	16 ng/mL
A/Shangdong/9/93	H3N2	16 ng/mL
A/Wisconsin/67/05	H3N2	31 ng/mL
A/Florida/04/06	Yamagata	8 ng/mL
B/Hong Kong/5/72	Victoria	31 ng/mL
B/Tokyo/53/99	Victoria	16 ng/mL
B/Victoria/504/00	Victoria	31 ng/mL

ANALYTICAL SPECIFICITY

The monoclonal antibodies used to detect SARS-CoV-2 specifically detect the virus' N protein.

The monoclonal antibodies used to detect Influenza A and B specifically detect a highly conserved epitope from the nucleoprotein of each virus.

DIAGNOSTIC SENSITIVITY AND SPECIFICITY

The Simple Ag SARS-CoV-2/Flu A+B test was evaluated by analysing different respiratory sample specimens.

The reference techniques were:

- RT-PCR CLART® PneumoVir 2(Genomica).
- PCR Cobas Liat Influenza A/B (Roche).
- Multiple Real Time PCR Kit for detection of 2019-CoV (XABT).
- Real SARS-CoV-2/Flu/RSV (Operon).

The results are shown below and separated by analyte.

A) Ag SARS-CoV-2 strip

The Ag SARS-CoV-2 test was evaluated both internally and externally at Hospital Clínico de Valencia (Spain).

A) Internal evaluation

390 respiratory samples were analysed, 117 of which were diagnosed as positive using the Multiple Real Time PCR Kit for detection of 2019-CoV (XABT) reference test.

The samples were classified as follows:

Samples	Overall	Oropharyngeal	Nasopharyngeal
Total	390	151	239
Total positives	117	66	51
Positives with Ct ≥ 30	31 (26.5 %)	19	12
Negatives	273	85	188

The following table contains the results.

593 samples		RT-PCR XABT		
Test	Result	Positive	Negative	Total
Simple Ag SARS-CoV-2/ Flu A+B	Positive	98	0	98
	Negative	19	273	292

The agreement of the Ag SARS-CoV-2 test results with RT-PCR were as follows:

Diagnostic specificity: **> 99.9 %**

Expressed as Ct values, the diagnostic sensitivity by viral load range is reported in the table below:

Specimen		Overall	Oropharyngeal	Nasopharyngeal
Diagnostic sensitivity %	Total	83.8	86.4	80.4
	Ct ≤ 30	90.6	89.1	92.7
	Ct ≤ 29	94.2	91.5	97.4
	Ct ≤ 28	95.2	93.2	97.4
	Ct ≤ 27	96.0	94.7	97.3

A) External evaluation

The external evaluation was performed at Hospital Clínico Universitario de Valencia (Spain). 167 nasopharyngeal respiratory samples from patients with symptoms compatible with COVID-19 were analysed, 53 of which were diagnosed as positive and 114 as negative, using the ThermoFisher Scientific RT-PCR TaqPath™ COVID-19 CE-IVD RT-PCR Kit reference test.

The following table contains the results.

167 samples		RT-PCR		
Test	Result	Positive	Negative	Total
Simple Ag SARS-CoV-2/ Flu A+B	Positive	32	1	33
	Negative	21	113	134

The agreement of the Ag SARS-CoV-2 test results with RT-PCR were as follows:

Diagnostic specificity: **> 99.1%**

Expressed as Ct values, the diagnostic sensitivity by viral load range is reported in the table below:

	Diagnostic sensitivity
All samples	60.4%
Ct ≤ 30	76.2%
Ct ≤ 29	80.0%
Ct ≤ 28	88.9%
Ct ≤ 27	88.9%
Ct ≤ 26	91.4%

B) Influenza strip

190 respiratory samples were analysed for Influenza A, 77 of which were diagnosed as positive. 145 samples were analysed for Influenza B, 14 of which were diagnosed as positive.

The results are shown below and separated by analyte, specimen, and reference technique.

Oropharyngeal swab samples in viral transport medium characterised by RT-PCR CLART® PneumoVir 2 (Genomica)

36 samples negative for Influenza A/B were analysed, along with 34 positive for Influenza A and 4 positive for Influenza B.

Influenza A		RT-PCR CLART® PneumoVir 2	
		Positive	Negative
Simple Ag SARS-CoV-2/ Flu A+B	Positive	20	0
	Negative	14	36
Sensitivity agreement = 58.8%			
Specificity agreement ≥ 99.9%			
PPV: 100 %			
NPV: 72.0%			

Influenza B		RT-PCR CLART® PneumoVir 2	
		Positive	Negative
Simple Ag SARS-CoV-2/ Flu A+B	Positive	2	0
	Negative	2	36
Sensitivity agreement = 50.0%			
Specificity agreement ≥ 99.9%			
PPV: 100%			
NPV: 94.7%			

Oropharyngeal swab samples in viral transport medium characterised by RT-PCR Real Time Cobas Liat Influenza A/B (Roche)

26 samples negative for Influenza A/B and 22 positive for Influenza A were analysed.

Influenza A		RT-PCR Cobas Liat Influenza A/B	
		Positive	Negative
Simple Ag SARS-CoV-2/ Flu A+B	Positive	14	1
	Negative	8	25
Sensitivity agreement = 63.6% Specificity agreement = 96.1% PPV: 93.6% NPV: 75.8%			

Influenza B		RT-PCR Cobas Liat Influenza A/B	
		Positive	Negative
Simple Ag SARS-CoV-2/ Flu A+B	Positive		0
	Negative		48
Specificity agreement ≥ 99.9% PPV: N/A NPV: 100%			

Oropharyngeal swab samples characterised by Real SARS-CoV-2/Flu/RSV (Operon)

47 samples negative for Influenza A/B were analysed.

Influenza A		RT-PCR	
		Positive	Negative
Simple Ag SARS-CoV-2/ Flu A+B	Positive		1
	Negative		46
Specificity agreement = 97.9% PPV: N/A NPV: 100%			

Influenza B		RT-PCR	
		Positive	Negative
Simple Ag SARS-CoV-2/ Flu A+B	Positive		0
	Negative		47
Specificity agreement ≥ 99.9% PPV: N/A NPV: 100%			

REPEATABILITY

Five replicates of a series of 1/2 dilutions of characterised internal standards were analysed (sensitivity curve), along with real samples established as NC ("negative sample"), LPC ("weak positive sample") and PC ("positive sample") for each analyte detected by the test, all measured on the same day by the same person. High repeatability was

obtained with differences between the replicates of the standard for each analyte equal to or below one 1/2 dilution and the results with the real samples were the same.

REPRODUCIBILITY

• **INTER-DAY PRECISION**

Over five consecutive days and in duplicate, the same Simple Ag SARS-CoV-2/Flu A+B test batch was measured for a sensitivity curve and a group of real samples established as NC ("negative sample"), LPC ("weak positive sample") and PC ("positive sample") for each of the analytes detected using the test. Sensitivity differences equal or less than one 1/2 dilution were obtained and the results with the real samples were the same.

• **INTER-OPERATOR PRECISION**

In triplicate, three persons measured a sensitivity curve and a group of real samples established as NC ("negative sample"), LPC ("weak positive sample") and PC ("positive sample") for each one of the analytes detected with the test. Sensitivity differences equal or less than one 1/2 dilution were obtained and the results with the real samples were the same.

• **INTER-BATCH PRECISION**

In parallel and in duplicate using three different test batches, a sensitivity curve was evaluated and a group of real samples established as NC ("negative sample"), LPC ("weak positive sample") and PC ("positive sample") for each of the analytes detected with the test. The analysis was performed by one person on the same day. Sensitivity differences equal or less than one 1/2 dilution were obtained and the results with the real samples were the same.

THE HOOK OR PROZONE EFFECT

The following was confirmed for the Simple Ag SARS-CoV-2/Flu A+B test:

1. The test does not present a Hook effect up to an N protein concentration of 800 ng/mL and up to an inactivated SARS-COV-2 (USA/WA1/2020) virus (ATCC NR-52287) concentration of 1×10^5 TCID₅₀/mL.
2. A Hook effect does not appear until an infectious viral load of 3.2×10^5 TCID₅₀/mL of the Influenza A/Virginia/ATCC1/2009 virus (ATCC Ref. VR-1736)
3. A Hook effect does not occur for Influenza B up to a concentration of 2.2×10^7 CEID₅₀/mL of the virus and Influenza B/Florida/04/06 (ATCC ref. VR-1804).

INTERFERING SUBSTANCES

The following substances showed no effect on the test results when added at the indicated concentrations:

Active ingredient	Concentration
Blood plasma	2%
Mucin	0.1%
Benzocaine	1.5 mg/mL
Hyaluronic acid and allantoin	1% v/v
Phenylephrine	15% v/v
Oxymetazoline	15% v/v
Momentasone	0.05 mg/mL
Fluticasone propionate	5% v/v
2,4-Dichlorobenzyl alcohol and 0.6 mg of amylmetacresol	5% m/v (50 mg/mL)
Sodium cromoglycate	15% v/v
Oseltamivir	5 mg/mL
Mupirocin	10 µg/mL
Tobramycin	4 µg/mL
Amoxicillin	1 mg/mL
Biotin	1 µg/mL
Alkalol (homeopathic)	1:10
Phenolic compounds	10% v/v
Ibuprofen	2.2 mg/mL
Paracetamol	1 mg/mL
Zanamivir	5 mg/mL
Ribavirin	1 mg/mL
Acetylsalicylic acid	0.3 mg/mL
Budesonide	10% v/v

CROSS-REACTIVITY

There is no evidence of cross-reactivity with the following micro-organisms:

Micro-organism	Analysed concentration
Bordetella pertussis	1 x 10 ⁶ CFU/mL
Haemophilus influenzae	1 x 10 ⁶ CFU/mL
Legionella pneumophila	1 x 10 ⁶ CFU/mL
Inactivated mycobacterium tuberculosis	1:10 dilution (0.27 g/mL)
Streptococcus pneumoniae	1 x 10 ⁶ CFU/mL
Streptococcus salivarius	1 x 10 ⁶ CFU/mL
Streptococcus pyogenes	1 x 10 ⁶ CFU/mL
Staphylococcus aureus	1 x 10 ⁶ CFU/mL
Staphylococcus epidermis	1 x 10 ⁶ CFU/mL
Candida albicans	1 x 10 ⁶ CFU/mL
Pseudomonas aeruginosa	1 x 10 ⁶ CFU/mL
Nasal wash	N/A
Inactivated coronavirus 229E	0.94 x 10 ⁴ TCID ₅₀ /mL
Inactivated coronavirus OC43	5 x 10 ⁴ TCID ₅₀ /mL

Viable coronavirus OC43	0.89 x 10 ⁴ TCID ₅₀ /mL
Inactivated coronavirus NL63	1.7 x 10 ⁴ TCID ₅₀ /mL
Inactivated MERS	1.7 x 10 ⁴ TCID ₅₀ /mL
Inactivated adenovirus 7A	1.5 x 10 ⁵ TCID ₅₀ /mL
Metapneumovirus TN/91-316	2.8 x 10 ⁵ TCID ₅₀ /mL
Parainfluenza 1 HPIV1/FRA/27344044/2007	0.89 x 10 ⁴ TCID ₅₀ /mL
Parainfluenza 2 (Greer)	0.89 x 10 ⁵ TCID ₅₀ /mL
Parainfluenza 3 (NIH 47885)	1.6 x 10 ⁶ TCID ₅₀ /mL
Parainfluenza 4B (19503)	5 x 10 ⁵ TCID ₅₀ /mL
Rhinovirus 48 (1505)	5 x 10 ⁵ TCID ₅₀ /mL
Enterovirus D 68 (USA/2018-23089)	2.8 x 10 ⁵ TCID ₅₀ /mL
RSV B	8.9 x 10 ⁴ TCID ₅₀ /mL

The Ag SARS-CoV-2 strip does not present cross-reactivity with Influenza A/B.

The Flu A+B strip does not present cross-reactivity with SARS-CoV-2.

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OPERON, S.A. - Camino del Plano, 19 -E-50410 Cuarte de Huerva - (Zaragoza) – SPAIN