



after bite
pain relief

BR 10 Insect bite healer



after bite
pain relief

For treating insect bites and stings



anti itching +
swelling

Can help combat itching and swelling



with
snap hook

Extra small design with snap hook
Always with you and ready to use - the companion
for everyday use, travel and sports



Heat can accelerate healing



Also suitable if you are pregnant as there are no chemical additives - only by means of warmth

Heat can accelerate healing

Quick-warming ceramic hot plate

3-second programm - also for sensitive skin

Medical device

Incl. 2 x 1.5V AAA batteries

Product measurements (LxWxH): 8.6 x 3.5 x 2.5 cm

Product weight: approx. 44.5 g (with batteries)

5-year warranty

Sales Unit: 12 / Shipping carton: 2 x 12

EAN no.: 4211125 603182

Item no.: 603.18