

BM 51
easyClip
beurer
medical

BM 51 easyClip - Upper arm blood pressure monitor


**german
engineering**



Innovative easyClip cuff: For upper arm circumferences from 22-42 cm



Inflation technology: Fast and convenient measurement, even during pumping-up



White illuminated, perfectly readable XL display



Fully automatic blood pressure and pulse measurement on the upper arm

Risk indicator: Classification of measurements with coloured scale

Arrhythmia detection: Warning in case of possible heart rhythm disturbances

Two user memories for 100 measured values each

Average value of all stored measurements and morning and evening blood pressure for the last 7 days

Can be used with mains part (available separately)

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Incl. 4 x 1.5V AA batteries and practical storage pouch

Medical device

Product measurements (LxWxH): 138 x 103 x 44 mm

Product weight: approx. 571 g (without batteries, with cuff)


Mains part available separately: Art. No. 071.95

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 4211125 65204 3

Item no.: 652.04



BM 51
easyClip
beurer
medical

easyClip – an innovative kind of cuff.

Elderly people in particular often have difficulties in putting the blood pressure cuff on the upper arm in the correct manner. Our innovative easyClip cuff can be put on with just one hand, which makes taking measurements even easier, quicker and more convenient.

clip & start



Clip the cuff on to your upper arm with just one hand...



...close the hook and loop fastener, ensuring a snug fit...



...and you're all set to measure your blood pressure on your upper arm.

easyClip cuff

Suitable for upper arm circumferences of 22-42 cm

