



BM 45

beurer
medical

BM 45 – Blood pressure monitor

 **german
engineering**

risk
indicator



Risk indicator:

Classification of measurements with coloured scale

xl display
+ light



White illuminated,

easy-to-read XL display



Fully automatic blood pressure and pulse measurement on the upper arm

Illuminated START/STOP button

Easy to use via sensor touch buttons

Arrhythmia detection: Warning in case of possible heart rhythm disturbance

Two user memories for 60 measured values each

With cuff for upper arm circumferences from 22-36 cm

Average value of all stored measurements and morning and evening blood pressure for the last 7 days

Can be used with mains part (available separately)

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Incl. 4 x 1.5V AA batteries and practical storage pouch

Medical device

Product measurements (LxWxH): 160 x 107 x 50 mm

Product weight: approx. 408 g (without batteries, with cuff)

XL cuff (35-44 cm) available separately: Art. No. 163.387

Mains part available separately: Art. No. 071.95

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 4211125 65806 9

Item no.: 658.06