



BM 28
Onpack

beurer
medical

BM 28 Onpack Blood pressure monitor

german
engineering



Resting indicator* for exact measurement results



With cuff for upper arm circumferences from 22-42 cm



Cuff position control**: For ensuring the cuff is correctly positioned on the upper arm



Fully automatic blood pressure and pulse measurement on the upper arm

Clear, easy-to-read display

Measurement reminder: 2 adjustable alarm times

Risk indicator: classification of measurements with coloured scale

Arrhythmia detection: Warning in case of possible heart rhythm disturbance

Four user memories for 30 measured values each

Average of all measurements

Date and time

Automatic switch-off

Incl. 4x1,5 V batteries AA and storage bag

Incl. mains adapter

Medical device

Mains operation:

Input: 100-240 V~; 50/60 Hz; max. 0.5 A

Output: 6V --- 600mA

Product measurements (LxWxH): 134 x 103 x 60 mm

Product weight: approx. 367 g (with cuff, without batteries)

5 year guarantee

Sales unit: 4 / Shipping carton: 6 x 4


EAN no.: 4211125 65830 4

Item no.: 658.30

****Cuff position control**



Tight
enough



Not tight
enough

658.30_BM28_Onpack-HSD_2019-08-20_01_DS_EN_BEU Subject to errors and changes

*Patent EP 1 673 009