



BC 44

beurer
medical


BC 44 wrist blood pressure monitor

 german
engineering



easy
to use

Easy to use via one single button



xl display
+ light

Blue illuminated, easy-to-read
XL display



Fully automatic blood pressure and pulse measurement on the wrist

With cuff for wrist circumferences from 12.5 - 21.5 cm

Arrhythmia detection: Warning in case of possible heart rhythm disturbances

Illuminated START/STOP button

Risk indicator: Classification of measurements with coloured scale

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Battery changing indicator

Incl. 2 x 1.5V AAA batteries and practical storage case

Medical device

Product measurements (LxWxH): 87.8 x 65.8 x 28.6 mm

Product weight: approx. 98 g (without batteries, with cuff)

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 4211125 65905 9

Item no.: 659.05