



FT 65 Multi-functional thermometer





Forehead-, ear- and surface temperature, high temperature alarm, date and time, 10 memory spaces



Measurement in seconds

6 in 1:

- Measurement within seconds in the ear
- Forehead temperature measurement in seconds
- Surface temperature of objects and liquids
- Optical high temperature alarm
- Date and Time
- 10 memory spaces

Suitable for the whole family

Large, easy-to-read display

Indication of normal or hightened temperature via coloured LED

- green = everything o.k.
- red = fever alarm (over 38° C / 100.4° F)

Measurement accuracy:

- Forehead: +/- 0.3 °C (+/- 0.5 °F)
- Ear: +/- 0.2 °C (+/- 0.4 °F)
- Object: +/- 1.5 °C (+/- 2.7 °F) with <30°C; +/- 5% with >30°C

Indication of °C and °F

Modern, handy design

Easy to clean

Automatic switch off

With battery changing indicator

Including 2 x 1.5V Battery AAA

Product measurements (L/W/H): $38.2 \times 138 \times 46.5 \text{ mm}$

Product weight: approx. 63g (without batteries)

Medical device

5 year guarantee

Sales unit: 4 / Shipping carton: 8 x 4

EAN-Nr.: 4211125 79514 6

Art.-Nr.: 795.14

195.14_FT65_2018-04-18_01_DS_EN_BEU Subject to errors and change