



ETM TESTMAGAZIN - URTEIL	
Beurer BM 54	
SEHR GUT	94,1 %
Blutdruckmessgerät im Online-Einzeltest www.etm-testmagazin.de	
	04/2019



BM 54

beurer
medical

BM 54 Bluetooth® - Upper arm blood pressure monitor

german
engineering



Universal cuff: also suitable for larger upper arm circumferences from 22-44 cm



Risk indicator: classification of measurements with coloured scale



- Fully automatic blood pressure and pulse measurement on the upper arm
- Innovative networking between smartphone and blood pressure monitor via Bluetooth® low energy technology
- Arrhythmia detection: Warning in case of possible heart rhythm disturbances
- Large, easy-to-read XL display
- Two user memories for 60 measured values each
- Average value of all stored measurements and morning and evening blood pressure for the last 7 days
- Date and time
- Automatic switch-off
- Signal in case of application error
- Automatic pressure preselection and deflation
- Incl. 4 x 1,5V AAA batteries and practical storage pouch
- Medical device
- Product measurements (LxWxH): 139 x 94 x 48 mm
- Product weight: approx. 364 g (without batteries, with cuff)
- 5 year guarantee
- Sales unit: 4 / Shipping carton: 4 x 4
- EAN no.: 4211125 65512 9
- Item no.: 655.12

655.12_BM54_2019-06-07_02_DS_EN_BEU Subject to errors and changes



ETW TESTMAGAZIN - URTEIL	
Beurer BM 54	
SEHR GUT	94,1 %
Blutdruckmessgerät im Online-Einzelttest www.etw-testmagazin.de	
	04/2019




BM 54

beurer
medical



CONNECTION



Innovative networking between smartphone and blood pressure monitor - keep track of your measurements at home and on the move. Easily synchronised with your smartphone via Bluetooth®.



Syncs with
• iOS ≥ 10.0
• Android™ ≥ 5.0
• Bluetooth® ≥ 4.0



List of supported devices

65512_BM54_2019-06-07_02_DS_EN_BEU Subject to errors and changes

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.

* tested online application "beurer HealthManager" Cloud
Can only be used with its full functionality with the „beurer HealthManager“ app.