





OT 20 Ovulation checking thermometer


 german
 engineering


Can be used with the "Ovy" app
 (manual recording of temperature)



Calculation of the cycle for
 natural family planning



Flexible measuring tip



Digital ovulation checking thermometer

Measurement of basal body temperature

Display of 2 digits after comma

For oral and vaginal measurement

Measurement accuracy:

± 0.05 °C: 35.00 °C - 38.00 °C at ambient temperature 18 °C - 28 °C;
 outside this range ± 0.1 °C

Measurement range: 32.00 - 42.99 °C

Waterproof

Can be disinfected

1 memory

Signal tone at end of measurement

Extra large display

With protection cap

Including 1x 1.5V Battery LR41

Battery changing indicator

Automatic switch off

Indication of °C

Product measurements (LxWxH): 142 x 23 x 13 mm

Product weight: approx. 12.9 g

Medical device

5 year guarantee

Sales unit: 20 / Shipping carton: 9 x 20

EAN no.: 4211125 79138 4

Item no.: 791.38