



BM 49

beurer
medical

BM 49 – Speaking upper arm blood pressure monitor

 **german
engineering**

speaking
function



Individually adjustable voice output in several languages (disengageable):

(a) DE, FR, IT, NL (c) RO, PL, CS, HU
(b) EN, ES, RU, GR (d) FA, AR, EN

xl display
+ light



Blue illuminated,
easy-to-read XL display

adjustable
volume



Adjustable volume



Fully automatic blood pressure and pulse measurement on the upper arm

Risk indicator: Classification of measurements with coloured scale

Arrhythmia detection: Warning in case of possible heart rhythm disturbance

Two user memories for 60 measured values each

With cuff for upper arm circumferences from 22–36 cm

Average value of all stored measurements and morning and evening blood pressure for the last 7 days

Can be used with mains part (available separately)

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Incl. 4 x 1,5V AA batteries and practical storage pouch

Medical device

Product measurements (LxWxH): 96 x 146 x 59 mm

Product weight: approx. 391 g (without batteries, with cuff)

XL cuff (35–44 cm) available separately: Art. No. 163.387

Mains part available separately: Art. No. 071.95

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: (a) 4211125 65628 7 (c) 4211125 65631 7

(b) 4211125 65640 9 (d) 4211125 65632 4

Item no.: (a) 656.28 (c) 656.31

(b) 656.40 (d) 656.32