



BM 96
CARDIO
beurer
medical

BM 96 **CARDIO** Blood pressure monitor with ECG function

 **german
engineering**

 **Bluetooth®**
beurer
CONNECT 



1-channel ECG device
for recording heart rhythm



Advanced arrhythmia detection: Detects atrial fibrillation
for stroke prevention



Resting indicator* for
exact measurement
results



Optimum control of the blood pressure and heart rhythm via app + web

Fully automatic blood pressure and pulse measurement on the upper arm

Risk Indicator: Classification of measurements

Easy-to-read XL display

Easy to use via 4 sensor touch buttons

Two user memories for 60 measured values each (With 30-second measurement. Deviation possible in the event of a longer measurement duration.)

Alarm function

Medical device

Date and time /automatic switch-off

Battery changing indicator

With cuff for upper arm circumferences from 22 - 42 cm

Incl. 4 x 1,5V AAA batteries and storage pouch

Product measurements (LxWxH): Blood pressure monitor: 120 x 144 x 53 mm
ECG stick: 120 x 22 mm

Product weight: Blood pressure monitor: ca. 490 g (without batteries, with cuff)
ECG stick: ca. 40 g

5 year guarantee

Sales unit: 4 / Shipping carton: 2 x 4

EAN no.: 4211125 65801 4

Item no.: 658.01

* Patent EP 1 673 009

658.01_BM96_2021-03-08_01_DS_EN_BEU Subject to errors and changes





BM 96
CARDIO


beurer
medical





CONNECTION













Innovative networking between smartphone and blood pressure monitor - keep track of your measurements at home and on the move. Easily synchronised with your smartphone via Bluetooth®.

works with





HealthManager Pro

works with




HealthManager






Synchs with



- iOS ≥ 12.0
- Android™ ≥ 8.0
- Bluetooth® ≥ 4.0
- ≥ Windows 8.1
- ≥ USB 2.0 (Type-A)



List of supported devices

658.01_BM96_2021-03-08_01_DS_EN_BEU Subject to errors and changes

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Android is a trademark of Google LLC.