




## OT 30 Ovulation checking thermometer


 german  
 engineering


Can be used with the "Ovy" app - Data transmission of temperature with *Bluetooth*® low energy technology



Calculation of the cycle for natural family planning



Documentation of more body signals



Digital ovulation checking thermometer

Measurement of basal body temperature

Display of 2 digits after comma

For oral and vaginal measurement

Illuminated LED ring signals *Bluetooth*® transmission

Measurement accuracy:

± 0.05 °C: 35.00 °C - 38.00 °C at ambient temperature from 18 °C to 28 °C;  
 outside this range ±0.1 °C

Measurement range: 32.00 - 42.99 °C

Flexible measurement tip

Waterproof

Can be disinfected

30 memory spaces

Signal tone at end of measurement (On/Off)

With protection cap

Including 1x Battery CR2032

Battery changing indicator

Automatic switch off

Indication of °C/°F

Product measurements (LxWxH): 105 x 33 x 17.6 mm

Product weight: approx. 20.1 g (without battery)

Medical device

5 year guarantee

Sales unit: 4 / Shipping carton: 12 x 4

EAN no.: 4211125 79106 3

Item no.: 791.06



german  
engineering



## CONNECTION



Calculation of the cycle for natural family planning together with the "Ovy" app - easy transmission of temperature with *Bluetooth*® low energy technology

works with



Syncs with

- iOS ≥ 8.0
- Android™ ≥ 4.4
- *Bluetooth*® 4.0

791.06\_OT30\_2018-07-02-03\_DS\_EN\_BEU Subject to errors and changes

The *Bluetooth*® word mark and logos are registered trademarks owned by *Bluetooth SIG, Inc.* and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.