



BM 95

beurer
medical

BM 95 Blood pressure monitor with ECG function

 german engineering

 **Bluetooth®**
beurer
CONNECT 



1-channel ECG device for recording heart rhythm



Advanced arrhythmia detection:
Detects atrial fibrillation and premature ventricular contractions



With cuff for upper arm circumferences from 22 - 42 cm



Optimum control of the blood pressure and heart rhythm via app + software

Fully automatic blood pressure and pulse measurement on the upper arm

Risk Indicator: Classification of measurements

Easy-to-read XL display

Two user memories for 60 measured values each

Alarm function

Medical device

Date and time /automatic switch-off

Battery changing indicator

Incl. 4 x 1,5V AAA batteries and storage pouch

Product measurements (LxWxH): Blood pressure monitor: 127 x 127 x 44 mm
ECG stick: 127 x 25 mm


Product weight: Blood pressure monitor: ca. 430 g (without batteries, with cuff)
ECG stick: ca. 40 g

5 year guarantee


Sales unit: 4 / Shipping carton: 2 x 4

EAN no.: 4211125 65834 2

Item no.: 658.34



CONNECTION



Can only be used with its full functionality with the „beurer CardioExpert“ app.
Free access to Beurer CardioExpert

658.34_BM95_2019-01-28_01_DS_EN_BEU Subject to errors and changes

The **Bluetooth®** word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Android is a trademark of Google LLC.