



BC 54

beurer
medical

BC 54 Bluetooth® Wrist blood pressure monitor

 german
engineering

 **Bluetooth®**
beurer
CONNECT 

inflation
technology



Fast and convenient measurement,
even during pumping up

risk
indicator




Classification of measurement
with coloured scale



Fully automatic blood pressure and pulse measurement on the wrist
With cuff for wrist circumferences from 13.5 - 21.5 cm
Clear, easy-to-read black display
Arrhythmia detection: Warning in case of possible heart rhythm disturbances
Two user memories for 60 measured values each
Average value of all stored measurements and morning and evening blood pressure for the last 7 days
Date and time
Automatic switch-off
Signal in case of application error
Automatic pressure preselection and deflation
Battery changing indicator
Incl. 2 x 1,5V AAA batteries and practical storage case
Medical device
Product measurements (LxWxH): 80.5 x 69.5 x 25 mm
Product weight: approx. 116 g (without batteries, with cuff)
5 year guarantee
Sales unit: 4 / Shipping carton: 4 x 4
EAN no.: 4211125 65054 4
Item no.: 650.54



CONNECTION



Innovative networking between smartphone and blood pressure monitor - keep track of your measurements at home and on the move. Easily synchronised with your smartphone via Bluetooth®.



Syncs with
• iOS ≥ 10.0
• Android™ ≥ 5.0
• Bluetooth® ≥ 4.0



List of supported devices

650.54_BC54_2019-09-19_01_DS_EN_BEU Subject to errors and changes

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.

* tested online application "beurer HealthManager" Cloud
Can only be used with its full functionality with the „beurer HealthManager“ app.