



german|engineering

## FT 100 Non-contact clinical thermometer

no  
contact

Contactless infrared clinical thermometer  
 Measures without skin contact - hygienic and safe

multi  
measuring

Measurement of forehead,  
 object and room temperature in seconds

no  
contact

Measurement distance sensor (signal tone/blue LED),  
 Measurement distance 4-6 cm



60 memory spaces (forehead measuring)

°C/°F can be set

Display illumination (On/Off)

Sound control (On/Off)

Automatic switch off

Incl. storage box

With battery changing indicator

Medical device

Incl. 2 x 1.5V batteries AAA

Product measurements (LxWxH): 199 x 41 x 68 mm

Product weight: approx. 82 g

3 year guarantee

Sales unit: 2 / Shipping carton: 8 x 2

EAN no.: 4211125 79505 4

Item no.: 795.05

≥ 38°C fever alarm    37,5-37,9°C elevated temperature    everything is OK

