

FT 09/1



FT 09/1 Thermometer white



beurer 370°



100% waterproof



Automatic switch off





With protective cover

Digital thermometer

Large Display

3-digit display

Audible signal when measuring stops

1 memory place

Measurement accuracy ± 0.1°C (between 35.5 and 42°C)

Measurement accuracy ± 0.2°C (between 32 - 35.4°C and 42.1 - 42.9°C)

100% waterproof

Can be disinfected

Without glass, mercury free

Incl. 1 x 1,5V battery LR41

Battery changing indicator

Product measurements: 129 x 20 x 17 mm

Product weight: approx. 10 g (including battery)

Medical device

5 year guarantee

Sales unit: 12 / Shipping carton: 12 x 12

EAN no.: 4211125 79115 5

Item no.: 791.15

791.15_FT09-1_white_2018-04-18_01_DS_EN_BEU Subject to errors and changes